

**Fax this form with the product order form to 877.342.8968**

Date: \_\_\_\_\_

Bill To (EasyStand Supplier Only)*
Purchase Order #:
Quote #:
Account #:
Ordered By:
Name of RTS:
Supplier Company:
Address:
City, State, Zip:
Country:
Phone:
Fax:
email:

Ship To
Name:
Address:
Address:
City, State, Zip:
Phone:
Tag For:

Parts Order
<input type="checkbox"/> Parts Order - Please refer to the off-unit order form
Serial Number _____ <i>(Serial number required for parts order)</i>
<i>Off Unit Price Lists &amp; Order Forms are available online at <a href="http://www.easystand.com/downloads">www.easystand.com/downloads</a></i>

Client Info
Height:
Weight:
Diagnosis:
Seat Depth (A):
Seat to Footplate (B):
Facility:
Therapist:

**Seat Depth (A):**  
Measure the distance between the back of the buttocks and the back of the knee.

**Seat to Footplate (B):**  
Measure the distance between the foot/shoe and the knee bend.



Notes:

\*Altimate Medical is the manufacturer of EasyStand products and does not sell directly to consumers. To buy an EasyStand, please contact a durable medical equipment supplier. If you have a supplier you are working with (for other medical equipment), find out if they are an EasyStand supplier. If not, use our supplier locator at [www.easystand.com](http://www.easystand.com) to find a qualified rehab equipment supplier in your area.

# EasyStand<sup>®</sup> Evolv

## Retail Price List & Order Form

Purchase Order #:

Account #:

Completed Bill To/Ship To Form Included

USA Effective January 1, 2012

**EasyStand<sup>®</sup> Evolv<sup>®</sup>**.....\$2,858

Choose a Size:


- PNG50162 Medium  
 PNG50084 Large

Choose actuator handle location:

- Right Side (standard)  
 Left Side

Approximate Height Range  
 Weight Limit  
 Seat Depth Range (seat pivot to hip pivot)  
 Seat to Footplate Range (from seat pivot)

	Medium	Large
Approximate Height Range	4'0"-5'6"	5'0"-6'2"
Weight Limit	200 lbs.	280 lbs.
Seat Depth Range (seat pivot to hip pivot)	14"-19"	16"-21"
Seat to Footplate Range (from seat pivot)	11"-18"	14"-21"



The Standard EasyStand Evolv is a sit to stand stander that includes:

- Two 5" Rear Locking Casters
- Two 5" Front Wheels
- Manual Hydraulic Actuator with Handle
- Planar Seat
- Flip-up Kneepad
- Multi-Adjustable Foot Plates
- Black Molded Tray with Chest Pad
- Transfer Handles
- Black Upholstery

### Frame Styles

- PNG50025 Mobile\*  
(self-propelled mobility, not available with independent kneepads).....\$1,266
- PNG50314 Front Swivel Casters†.....\$153
- PNG50417 Swing Away Front†.....\$231
- PNG50428 Swing Away Front for Shadow Tray†.....\$180

### Footplates and Seats

- PNG30030 Foot Straps.....\$ 37
- PNG30031 Secure Foot Straps.....\$ 99
- PNG30269 Transfer Seat (Medium only).....\$118
- PNG30270 Transfer Seat (Large only).....\$118
- PNG50324 Rotating Seat (Medium only).....\$206
- PNG50318 Rotating Seat (Large only).....\$206
- PNG30183 Seat with Roho<sup>®</sup> Insert (Medium only).....\$524
- PNG30023 Seat with Roho<sup>®</sup> Insert (Large only).....\$524
- P80246 Seat Angle Locator.....\$ 30
- PNG50161 Hip Supports-Medium.....\$312
- PNG50198 Hip Supports-Large  
(limited adjustment with Shadow Tray).....\$312
- PNG50384 Independent Kneepads†.....\$418
- PNG50385 Independent Roho<sup>®</sup> Kneepads†.....\$524

### Tray/Front Supports (back support required for all Shadow Trays)

- PNG50471 Shadow Tray 7"†  
(7.25" H from seat, 3.5"-15" D from back).....\$524
- PNG50474 Deep Shadow Tray 7"†  
(7.25" H from seat, 7.5"-19" D from back).....\$524
- PNG50130 Shadow Tray 10"†  
(10" H from seat, 3.5"-15" D from back).....\$524
- PNG50334 Deep Shadow Tray 10"†  
(10" H from seat, 7.5"-19" D from back).....\$524
- PNG50026 Shadow Tray 12"†  
(12.25" H from seat, 3.5"-15" D from back, Large only).....\$524
- PNG50332 Deep Shadow Tray 12"†  
(12.25" H from seat, 7.5"-19" D from back, Large only).....\$524
- PNG50248 Clear Tray Top.....\$139
- PNG50194 Oversized Clear Tray†  
(not available with Shadow Tray).....\$237
- PNG50368 Black Molded Angle Adjustable Tray\*.....\$301
- PNG50366 Clear Angle Adjustable Tray\*.....\$437
- PNG50327 Clear Angle Adjustable Tray for Mobile\*.....\$301
- PNG50295 No-Tray Chestpad\*.....\$217
- PNG50042 Hand Grips  
(clear tray required).....\$131
- P5030 Arm Troughs †  
(clear tray required).....\$275
- PNG30000 Large Contoured Chest Pad  
(standard with shadow tray).....\$121
- PNG30295 Padded Tray Cover†  
(not available with oversized clear tray).....\$194

\* not available with Shadow Tray or Swing Away Front  
 † not available with Mobile

### Lift Mechanisms

- PNG50010 Pow'r Up Lift  
(replaces hydraulic lift, battery powered).....\$1,054
- PNG50328 Adjustable/Removable Actuator Handle.....\$195
- PNG50044 Quad Grip Handle Extension  
(not available with Pow'r Up).....\$131
- PNG50070 T-Style Handle Extension  
(not available with Pow'r Up).....\$131

### Back / Upper Body Supports

- PNG50066 Flat Back 13".....\$572
- PNG50067 Flat Back 19".....\$604
- PNG50062 Contoured Back 13".....\$630
- PNG50068 Contoured Back 19".....\$662
- PNG50171 Removable Flat Back 13".....\$657
- PNG50173 Removable Flat Back 19".....\$689
- PNG50172 Removable Contoured Back 13".....\$715
- PNG50174 Removable Contoured Back 19".....\$747

**ATTENTION** - All options below require one of the back supports above

- PNG30029 Velcro<sup>®</sup> Positioning Belt.....\$ 47
- PNG30028 Positioning Belt with Airline Style Buckle.....\$131
- PNG50192 Chest Strap-Small.....\$121
- PNG50193 Chest Strap-Large.....\$121
- P600002 Back Angle Adjustment Wheel.....\$ 89

**ATTENTION** - All options below require the Accessories Mounting Bracket

- PNG50034 Accessories Mounting Bracket.....\$ 57
- PNG30170 X-Style Chest Vest-Medium.....\$237
- PNG30056 X-Style Chest Vest-Large.....\$237
- PNG30171 Y-Style Chest Vest-Medium.....\$237
- PNG30021 Y-Style Chest Vest-Large.....\$237
- PNG11212 High Mount Chest Vest Bracket.....\$ 66
- PNG50224 Lateral Supports 9"-15" W.....\$312
- PNG50176 Lateral Supports 12"-18" W.....\$312
- PNG50240 Head Support 21"-30" H.....\$312
- PNG50037 Head Support 24"-33" H.....\$312
- PNG50040 Push Handles.....\$131

**All options ordered on-unit are in place of similar standard items.**

### Option Reminder

- ✓ Back Option      ✓ Tray Option      ✓ Accessories Mounting Bracket

**Fax this form to 877.342.8968**

**AltimateMedical**

Ph: 800.342.8968 or 507.697.6393  
 Fax: 877.342.8968 or 507.697.6900  
 www.easystand.com