

A New Approach to Standing Kids



The EasyStand Bantam is new type of stander for kids! It is the only standing frame to combine a sit-to-stand transition with the added option of a supine stander.

The Bantam stander makes it safer and easier to transfer and stand your child, so your child will stand more often.

By alternating a variety of positions, your child can tolerate standing for a longer period of time, maximizing the health benefits.



Lily Harris  
Cerebral Palsy  
40" tall (102cm)

## Bantam Specifications

|                                 | Extra Small             | Small                  |
|---------------------------------|-------------------------|------------------------|
| <b>Height Range</b>             | 28" - 40"<br>71-102cm   | 36" - 54"<br>91-137cm  |
| <b>Weight Limit</b>             | 50 lbs.<br>23kg         | 100 lbs.<br>45kg       |
| <b>Seat Depth Range</b>         | 7" - 12"                | 11" - 16"              |
| -seat pivot to back surface     | 18-31cm                 | 28-41cm                |
| <b>Seat/Transfer Height</b>     | 19"                     | 19"                    |
| -from floor                     | 48cm                    | 48cm                   |
| <b>Seat to Foot Plate Range</b> | 4.5" - 16.5"            | 4.5" - 16.5"           |
| -from front of seat             | 11-42cm                 | 11-42cm                |
| <b>Kneepad Depth Range</b>      | 2.5" - 8.5"             | 2.5" - 8.5"            |
| -from front of seat             | 6-22cm                  | 6-22cm                 |
| <b>Knee Width</b>               | 5.5"<br>14cm            | 5.5"<br>14cm           |
| <b>Back Height Range</b>        | 11" - 13"               | 12" - 14"              |
| -from seat                      | 28-33cm                 | 30-36cm                |
| <b>Weight of Basic Unit</b>     | 54 lbs.<br>24.5kg       | 54 lbs.<br>24.5kg      |
| <b>Frame Footprint</b>          | 23" x 33.5"<br>58x85cm  | 23" x 33.5"<br>58x85cm |
| <b>Upholstery</b>               | Black Nylon/Neoprene    |                        |
| <b>Frame Color</b>              | Red or Blue Powder Coat |                        |
| <b>Lifting Mechanism</b>        | Gas Spring Lift         |                        |

*"Great for school, the ease of use will promote staff to utilize equipment with Emma. Provides weight bearing on hips & legs. Transition from sitting to standing is a good sensory experience and standing at eye level is great for my daughter's vision."*  
**Parent of Emma G. - Phoenix, AZ, USA**

## Unique Bantam Options

The Bantam meets the needs of kids of various abilities, with accessories like the self-propelled Mobile, Shadow Tray, supine positioning and numerous support options.

To gain the greatest benefits, a child with a disability should start standing around the same time as their able-bodied peers, which is usually around one year of age.

The EasyStand Bantam allows a child to start their standing program from supine, then progress to the preferred sit-to-stand method of standing as they grow heavier or gain more control. The EasyStand Bantam is available with a growth kit to fit kids from early intervention through elementary school.



Arianna Saldana  
Cerebral Palsy  
36" tall (92cm)

Bantam Mobile



## Sit. Stand. Smile. Repeat.

**Transfers** - Your child can be transferred into the Bantam stander with a lateral seat-to-seat transfer, supine flat-to-load transfer, or supine flat-to-load with flexion transfer.

**Transitions** - The natural sit-to-stand transition helps minimize your child's anxiety. A caregiver helps your child into a standing position with a gas spring lift assist, or the child can achieve independent standing with the manual hydraulic pump (optional) or Pow'r Up Lift (optional).

**Positions** - By combining the sit-to-stand function of the Bantam stander with the supine option, positions from sitting to standing, to supine and everything in between are achievable. Even children who have contractures can bear weight, and often progress to a more upright standing position.



90/90 Sitting



Anterior Tilt Sitting



Upright Standing



Supine with Hip-Knee Flexion



Supine Standing



Supine Lying

# Standard Features



## EasyStand Bantam



### #PT50001 EasyStand Bantam - Extra Small (XS)

### #PT50002 EasyStand Bantam - Small (S)

The EasyStand Bantam is a sit-to-stand stander for kids. It includes foot plates, kneepads, gas spring assist lift with foot pedal, planar seat and back, black upholstery, two front wheels, two rear locking casters and chest strap. It is available with a red or blue frame (tray not included).

# Options

## EasyStand Bantam



### #PT50009 Supine Positioning

Allows the Bantam to be used as a supine stander. Child can be transferred into stander in flat-to-load or seated position. Adjusts in 0-90° supine. Head support required.



### #PT50150 Mobile

Allows child to self-propel while sitting or standing. Push rims move forward for transferring and can be re-positioned to child's reach. Mobile footprint 24"x37" (61x94cm). Not available with Shadow Tray, Supine Positioning, or Front Swivel Casters.



### #PT50116 Front Swivel Casters

Locking casters allow for easier steering and movement of stander from front or back of unit. Adds 5" width to the front of unit. Not available with Mobile option.



### #PT50022 Pow'r Up Lift

Brings the child to the standing position with a touch of a button. Includes a light touch controller, rechargeable battery, Linak® actuator, and an emergency release.



### #PT50020 Hydraulic Pump w/Removable Pump Handle

Allows the child to independently transition to standing with a manual hydraulic pump. Handle is adjustable and removable.



### #PNG50044 Quad Grip Handle Extension

Attaches to the pump handle, for kids with limited dexterity. Adjusts for various hand sizes. Hydraulic pump with removable handle required.



### #PNG50070 T-Style Handle Extension

Attaches to the pump handle, for kids with limited range of motion. Makes the first few pumps easier and releasing more accessible. Adds additional 6" (15cm) in each direction. Hydraulic pump with removable handle required.



### #PT50140 Assist Handle Extension

Assists caregiver in moving seat from sit-to-stand, and stand to sit, by providing an extension to the current handle. Useful with the gas spring assist lift when the child is very lightweight.



### #PT50080 Foot Straps

D-Ring with Velcro adjustment holds feet in alignment with multiple attachment slots for proper foot positioning.



### #PT50082 Secure Foot Straps

D-Ring with Velcro adjustment holds feet securely in place at ankles and toes with multiple attachment slots for proper foot positioning.



### #PT50132 Multi-Adjustable Footplates (XS)

Deluxe footplates are adjustable in forward/aft +/- 3/4" (2cm), plantar/dorsi +/- 20°, toe-in/toe-out, swivel 30° left/right, in addition to height. Extra small size is 3 3/8" (9cm) wide by 7 7/8" (20cm) long.



### #PT50006 Multi-Adjustable Footplates (S)

Deluxe footplates are adjustable in forward/aft +/- 3/4" (2cm), plantar/dorsi +/- 20°, toe-in/toe-out, swivel 30° left/right, in addition to height. Small size is 4 7/8" (10cm) wide by 9 3/4" (12cm) long.



### #PT50094 Multi-Adjustable Kneepads (XS)

Deluxe kneepads are adjustable in swivel, height +/- 1" (2.5cm), width 5 1/2"-10" (14-25cm), depth range of 2 1/2"-8 1/2" (6-22cm) Inside access of pad is 2 1/2" (6cm) wide by 5" (13cm) height.



### #PT50004 Multi-Adjustable Kneepads (S)

Deluxe kneepads are adjustable in swivel, height +/- 1" (3cm), width 5 1/2"-10" (14-25cm), depth range of 2 1/2"-8 1/2" (6-22cm) Inside access of pad is 3 1/2" (9cm) wide by 5" (13cm) height.



**#PT50074 Black Molded Shadow Tray**  
Shadow tray provides anterior support as child moves from sitting to standing. Adjustable in depth, height and tray angle. Tray size is 19"x24" (48x61cm).



**#PT50076 Clear Shadow Tray**  
Shadow tray provides anterior support as child moves from sitting to standing. Adjustable in depth, height and tray angle. Tray size is 19"x24" (48x61cm).



**#PNG50042 Hand Grips**  
Provide additional forearm support and alignment. Attach to clear acrylic tray with suction cups.



**#PT0156 Elbow Pad-9.25" Cutout**  
Elbow/Forearm pad provides additional padded support surface on the black molded or clear tray. 9¼"(24cm) cutout width. Overall width 21½"(55cm), pad extends back 3¾"(10cm) posterior from current table edge for additional upper body/extremity support.



**#PT50160 Elbow Pad-11" Cutout**  
Elbow/Forearm pad provides additional padded support surface on the black molded or clear tray. 11"(28cm) cutout width. Overall width 21½"(55cm), pad extends back 3¾"(10cm) posterior from current table edge for additional upper body/extremity support.



**#PT50044 Contoured Seat (XS)**  
**#PT50048 Contoured Seat (S)**  
Contoured positioning seat features form-to-fit upholstery with built in hip support.



**#PT30058 Easy-Wash Seat Cover/Planar (XS)**  
**#PT30060 Easy-Wash Seat Cover/Planar (S)**  
**#PT30059 Easy-Wash Seat Cover/Contoured (XS)**  
**#PT30061 Easy-Wash Seat Cover/Contoured (S)**  
Tek Stretch seat cover can be easily wiped clean or laundered. Protects padded seat from spills or accidents.



**#PT50052 Contoured Back (XS)**  
**#PT50056 Contoured Back (S)**  
Contoured positioning back features form-to-fit upholstery with built in lateral support.



**#PT50098 Hip Supports (XS)**  
Support the hips in symmetrical alignment from sitting to standing. Adjustable from 6"- 9" (15-23 cm) width with a 2" (5cm) height range. Pad size is 3½"x 3½" (9x9cm). Not available with contoured seat.



**#PT50007 Hip Supports (S)**  
Support the hips in symmetrical alignment from sitting to standing. Adjustable from 8"-12" (20-30cm) width with a 2" (5cm) height range. Pad size is 3½"x 3½" (9x9cm). Not available with contoured seat.



**#PT50008 Lateral Supports**  
Provides additional midline positioning for children with limited trunk control. Width adjustment range is 8"-12" (20-30cm). Pad size is 4"W x 6"H (10x15cm). Not available with contoured back.



**#PT50124 Positioning Belt (XS)**  
**#PT50090 Positioning Belt (S)**  
2" wide (5cm) positioning belt has a D-ring velcro closure, provides hip stability for child. Fits hip circumference up to 23" (58cm) for extra small size and up to 29" (74cm) for small size.



**#PT50128 High Mount Chest Vest Bracket**  
Moves mounting points for chest vest higher to secure shoulders in position. Chest vest not included.



**#PT50086 Chest Vest (XS)**  
**#PT50088 Chest Vest (S)**  
Neoprene with camlock adjustment provides anterior trunk support for the child. Extra small size is 9½"L x 9"W (24x23cm). Small size is 11"L x 9½"W (28x24cm).



**#PT50062 Head Support (XS)**  
**#PT50064 Head Support (S)**  
Contoured head support helps to maintain optimal positioning of the head in midline. Height range from seat is 15"-20" (38-51cm), height range with inverted bracket is 17"-22" (43-56cm). Depth range is -½" to -3" (-1 to -8cm) from backrest.



**#PNG50040 Push Handles**  
Push handles make moving the stander from place to place easier for the caregiver.



**#T50118 Growth Kit - Planar Seating**  
Converts a Bantam Extra Small to a Bantam Small to accommodate for a growing child. Includes a planar seat, planar back and kneepads.



**#T50120 Growth Kit - Contoured Seating**  
Converts a Bantam Extra Small to a Bantam Small to accommodate for a growing child. Includes contoured seat, contoured back and kneepads.



**#P80246 Seat Angle Locator**  
Helps therapists monitor standing progress. Attaches under seat magnetically.



**#P80809 Tool Pouch**  
Holds the assembly tools provided with the stander. Attaches to unit with velcro straps. Tools not included.