

Essential Features, Economical Price, Easily Attainable

In response to supplier demand, the Evolv E3 base price is low, so you can “build” on the options your clients need while staying within funding guidelines.

The EasyStand Evolv E3 is only \$2,195 compared to \$2,775 for the Evolv. To lower the base price of the Evolv E3, standard features focus on the basic essentials.

Features of the Evolv E3 are different than the Evolv, including a no-tray chest pad that replaces the black molded tray and footplates that are height adjustable only.

The Evolv E3 is available in Youth or Adult sizes and offers all of the options and modularity of the Evolv. Please refer to the EasyStand Evolv E3 pricelists for details.

EasyStand Evolv E3
MSRP \$2,195

The footplates are adjustable in height only.

The black molded tray is replaced by a no-tray chest pad.

Casters are included!

EasyStand Evolv E3 - \$2,195

EasyStand standing frames are typically billed under the HCPCS payer code **E0637** (combination sit to stand system, any size including pediatric, with seat lift feature, with or without wheels) or **E1399** (miscellaneous) in most states.

EasyStand® Evolv E3

Retail Price List & Order Form

USA Effective January 1, 2011

Purchase Order #:

Account #:

Completed Bill To/Ship To Form Included

EasyStand® Evolv® E3.....\$2,195

The EasyStand Evolv E3 is different from the EasyStand Evolv. The standard features on the E3 model differ from those offered on the EasyStand Evolv. Before ordering the E3 model be certain that you know the differences and are ordering the product that best fits your client's needs.


Choose a Size:

- PNG50410 Medium (formerly Youth)
- PNG50400 Large (formerly Adult)

Choose actuator handle location:

- Right Side (standard)
- Left Side

	Medium	Large
Approximate Height Range	4'0"-5'6"	5'0"-6'2"
Weight Limit	200 lbs.	280 lbs.
Seat Depth Range (seat pivot to hip pivot)	14"-19"	16"-21"
Seat to Footplate Range (from seat pivot)	11"-18"	14"-21"



The Standard EasyStand Evolv E3 is a sit to stand stander that includes:

- Two 5" Rear Locking Casters
- Two 5" Front Wheels
- Manual Hydraulic Actuator with Handle
- Planar Seat
- Flip-up Kneepad
- Height-only Adjustable Foot Plates
- No-tray Chest Pad
- Transfer Handles
- Black Upholstery

Frame Styles

- PNG50025 Mobile* (self-propelled mobility, not available with independent kneepads).....\$1,230
- PNG50314 Front Swivel Casters†.....\$149
- PNG50417 Swing Away Front† (tray required).....\$225
- PNG50428 Swing Away Front for Shadow Tray† (Shadow tray required).....\$175

Footplates and Seats

- PNG50414 Multi-Adjustable Footplates (Medium only).....\$180
- PNG50411 Multi-Adjustable Footplates (Large only).....\$180
- PNG30030 Foot Straps.....\$ 36
- PNG30031 Secure Foot Straps.....\$ 97
- PNG30269 Transfer Seat (Medium only).....\$115
- PNG30270 Transfer Seat (Large only).....\$115
- PNG50324 Rotating Seat (Medium only).....\$200
- PNG50318 Rotating Seat (Large only).....\$200
- PNG30183 Seat with Roho® Insert (Medium only).....\$509
- PNG30023 Seat with Roho® Insert (Large only).....\$509
- P80246 Seat Angle Locator.....\$ 30
- PNG50161 Hip Supports-Medium.....\$303
- PNG50198 Hip Supports-Large (limited adjustment with Shadow Tray).....\$303
- PNG50384 Independent Kneepads†.....\$406
- PNG50385 Independent Roho® Kneepads†.....\$509

Front Mounted Trays (not available with Shadow Tray)

- PNG50393 Black Molded Tray†.....\$489
- PNG50394 Clear Tray†.....\$630
- PNG50395 Clear Tray for Mobile.....\$630
- PNG50396 Oversized Clear Tray†.....\$746
- PNG50397 Black Molded Angle Adjustable Tray*†.....\$782
- PNG50398 Clear Angle Adjustable Tray*†.....\$915
- PNG50399 Clear Angle Adjustable Tray for Mobile*.....\$915

Shadow Trays (back support required for all Shadow Trays)

- PNG50476 Shadow Tray 7"† (7.25" H from seat, 3.5"-15" D from back).....\$1,024
- PNG50478 Deep Shadow Tray 7"† (7.25" H from seat, 7.5"-19" D from back).....\$1,024
- PNG50401 Shadow Tray 10"† (10" H from seat, 3.5"-15" D from back).....\$1,024
- PNG50402 Deep Shadow Tray 10"† (10" H from seat, 7.5"-19" D from back).....\$1,024
- PNG50403 Shadow Tray 12"† (12.25" H from seat, 3.5"-15" D from back).....\$1,024
- PNG50404 Deep Shadow Tray 12"† (12.25" H from seat, 7.5"-19" D from back).....\$1,024

Tray Accessories (tray required)

- PNG50347 Clear Top for Shadow Tray (must order Shadow Tray).....\$135
- PNG50042 Hand Grips (clear tray required).....\$128
- P5030 Arm Troughs† (clear tray required).....\$267
- PNG30000 Large Contoured Chest Pad (standard with shadow tray).....\$118

* not available with Shadow Tray † not available with Mobile

Lift Mechanisms

- PNG50010 Pow'r Up Lift (replaces hydraulic lift, battery powered).....\$1,024
- PNG50328 Adjustable/Removable Actuator Handle.....\$190
- PNG50044 Quad Grip Handle Extension (not available with Pow'r Up).....\$128
- PNG50070 T-Style Handle Extension (not available with Pow'r Up).....\$128

Back/Upper Body Supports

- PNG50066 Flat Back 13".....\$556
- PNG50067 Flat Back 19".....\$587
- PNG50062 Contoured Back 13".....\$612
- PNG50068 Contoured Back 19".....\$643
- PNG50171 Removable Flat Back 13".....\$638
- PNG50173 Removable Flat Back 19".....\$669
- PNG50172 Removable Contoured Back 13".....\$695
- PNG50174 Removable Contoured Back 19".....\$726

ATTENTION- All options below require one of the back supports above

- PNG30029 Velcro® Positioning Belt.....\$ 46
- PNG30028 Positioning Belt with Airline Style Buckle.....\$128
- PNG50192 Chest Strap-Small.....\$118
- PNG50193 Chest Strap-Large.....\$118
- P600002 Back Angle Adjustment Wheel.....\$ 87

ATTENTION- All options below require the Accessories Mounting Bracket

- PNG50034 Accessories Mounting Bracket.....\$ 56
- PNG30170 X-Style Chest Vest-Medium.....\$231
- PNG30056 X-Style Chest Vest-Large.....\$231
- PNG30171 Y-Style Chest Vest-Medium.....\$231
- PNG30021 Y-Style Chest Vest-Large.....\$231
- PNG11212 High Mount Chest Vest Bracket.....\$ 65
- PNG50224 Lateral Supports 9"-15" W.....\$303
- PNG50176 Lateral Supports 12"-18" W.....\$303
- PNG50240 Head Support 21"-30" H.....\$303
- PNG50037 Head Support 24"-33" H.....\$303
- PNG50040 Push Handles.....\$128

All options ordered on-unit are in place of similar standard items.

Option Reminder

- ✓ Tray
- ✓ Back
- ✓ Multi-Adjustable Footplates/Kneepads
- ✓ Accessories Mounting Bracket

Fax this form to 877.342.8968

AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

www.easystand.com

Standing technology should only be used under the guidance of a physician with recommendations for standing program protocol and any medical precautions. Standing programs should be monitored by the attending therapist. AMI maintains a policy of continual product improvement and reserves the right to change features, specifications, and prices without prior notification. Check with AMI for latest information. FORM ESPL-E3-0111 Copyright © 2011 Altimate Medical, Inc. All rights reserved. U.S. & foreign patents pending. Printed in the U.S.A.

Fax this form with the product order form to 877.342.8968

Date: _____

Bill To (EasyStand Supplier Only)

Purchase Order #:

Quote #:

Account #:

Ordered By:

Name of RTS:

Supplier Company:

Address:

City, State, Zip:

Country:

Phone:

Fax:

email:

Ship To

Name:

Address:

Address:

City, State, Zip:

Phone:

Tag For:

Parts Order

Parts Order - Please refer to the off-unit order form

Serial Number _____
(Serial number required for parts order)

Off Unit Price Lists & Order Forms are available online at www.easystand.com/downloads

Client Info

Height:

Weight:

Diagnosis:

Seat Depth (A):

Seat to Footplate (B):

Facility:

Therapist:

Seat Depth (A):

Measure the distance between the back of the buttocks and the back of the knee.

Seat to Footplate (B):

Measure the distance between the foot/shoe and the knee bend.



Notes:

AltimateMedical

PO Box 180, 262 West 1st Street, Morton, MN 56270 USA

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

info@easystand.com www.easystand.com