

EasyStand® StrapStand



Wheel Up and Stand



The EasyStand StrapStand is a sit-to-stand stander that fits the needs of wheelchair users of many diagnoses, from low level spinal cord injuries to multiple sclerosis. It is often used in home standing programs, as well as multi-user settings, such as nursing homes or assisted living centers.

Newly redesigned in 2011, the EasyStand StrapStand is now available with more options than ever, many that are similar to the EasyStand Evolv.

The most popular option on the StrapStand continues to be the Swing-Out Legs, which are now easier to operate. With the simple shift of the lever, the caregiver can swing the legs of the stander in/out, similar to how they would operate legs on a patient lift.



Melanie Peck-Mosen
Multiple Sclerosis
5'2" (157cm)

"Outstanding! I love the EasyStand StrapStand because my wheelchair can be brought flush to the stander and it was cost effective. It promotes strength in my lower extremities, builds tolerance, and increases balance for proper postural control & visual attention to a task."
Sandy O. - Elkhorn, WI

"The StrapStand saves the backs of employees assisting patients to the standing position. We also like the portability because we are able to take it to patients' rooms. We use it daily with a variety of residents."
Kim B. - Wyandotte, MI, USA



Stand From Your Wheelchair

No Transfer Needed - If you are not capable of transferring, or patient lifting is a concern for staff or caregivers, the StrapStand may fit your needs. The StrapStand lifts you directly from your wheelchair, bed, or other seated surface into the standing position.

Transitioning - To move from your wheelchair to standing, position the lifting strap behind your hips and pull back on the handle to operate the manual hydraulic pump. The difference between the transition of the StrapStand and other strap-style standers is how it lifts, which affects the comfort of the knees. With the EasyStand StrapStand, the lifting arms gently lift you up to standing (rather than pushing you into the kneepad) which greatly reduce the pressure on your knees.

Standing - The StrapStand can provide the maximal over-center stretch needed for the optimum standing posture.



Carla Cardinal
Multiple Sclerosis
5'0" (152cm)

StrapStand Specifications

Height Range	5' - 6'5"	152-195cm
Weight Limit	350 lbs.	159kg
Kneepad Height Range		
-from foot plate	18" -21"	45-53cm
Kneepad Width		
-center to center	9"	23cm
Kneepad Depth Range		
-from back of heel	6" - 9.5"	15-24cm
Kneepad Size	9" x 15.5"	23x39cm
Front Pad Depth Range		
-from back of heel	2" -12.5"	5-32cm
Table Height Range		
-from foot plate	40" -56"	102-142cm
Inside Width / Wheelchair Access	22.75"	58cm
Footprint	31.5" x 27.5"	80x70cm
Frame Color	Charcoal	
Upholstery Color	Graphite Gray	
Weight of Basic Unit	135 lbs.	61kg

Standard Features

EasyStand StrapStand



#2100 EasyStand StrapStand

The EasyStand StrapStand is a sit-to-stand stander that includes two 3" locking casters, two 4" wheels, foot plates, kneepad, no-tray chest pad, manual hydraulic actuator with handle, black upholstery and charcoal frame (lifting strap not included).



#80811 Lifting Strap - Small
#80719 Lifting Strap - Large
#80962 Lifting Strap - XL

Padded lifting strap is available in three sizes. Small strap is 10" x 26" (25x66cm). Large strap is 12" x 29" (30x74cm). XL strap is 14" x 32" (36x81cm).



#80882 Adjustable Lifting Strap - Small
#80881 Adjustable Lifting Strap - Large
#80963 Adjustable Lifting Strap - XL

Adjustable lifting strap is useful for individuals using fixed arm wheelchairs. Small strap is 10" x 26" (25x66cm). Large strap is 12" x 29" (30x74cm). XL strap is 14" x 32" (36x81cm). Range on all three sizes is 0" - 8" (0-20cm).



PNG50393 Black Molded Tray

Black molded tray provides anterior support and work surface. Adjustable in height and depth. Tray size is 21" x 24" (53x56cm).



PNG50394 Clear Tray

Clear Acrylic tray provides anterior support and work surface. Adjustable in height and depth. Tray size is 21" x 24" (53x61cm).



PNG50398 Clear Angle Adjustable Tray
PNG50397 Black Molded Angle Adjustable Tray

Adjusts from 0°-35° to provide forearm support and angle activities closer to the user.



PNG50042 Hand Grips

Provide additional arm and hand support and alignment. Easily movable for precise arm positioning. Attaches to clear acrylic tray with suction cups.



P5030 Arm Troughs

Provide additional forearm support and alignment. Easily movable for precise arm positioning. Attach to clear acrylic tray with suction cups.



PNG30000 Large Contoured Chest Pad

Attaches to tray or no-tray bracket to provide anterior support with a minimal lateral hug. Replaces standard chest pad. Oval pad is 9½" x 15" (24x38cm).



P82090 Multi-Adjustable Footplates

Adjust independently in height from 1" - 5" off floor (2.5-13cm), -20°/+20° plantar/dorsi flexion adjustment, and +15° toe out adjustment.



P82144 XL Footplate

Accommodates larger feet. Each foot placement area is 7" x 14" (18x36cm). Not available with secure foot straps.



PNG30030 Foot Straps

Velcro® and d-ring adjustment hold feet into alignment with multiple attachment slots on footplates for proper foot positioning. Straps fit foot height of 4" - 6½" (10-17cm).



#PNG30031 Secure Foot Straps

Velcro® and D-ring adjustment holds feet securely in place at ankles and toes. Straps fit foot height of 4" - 6½" (10-17cm).



P82101 Independent Kneepads

Accommodates hip or knee contractures and leg length discrepancies. Kneepad width range from center to center is 8"-12" (20-30cm). Kneepad depth range from back of heel is 3½" - 8" (9-20cm). Kneepad height range from footplate to top of pad is 17" - 22½" (43-57cm).



P82189 Independent Roho® Kneepads

Accommodates hip or knee contractures and leg length discrepancies with Roho® Dry Floatation® technology. Kneepad width range from center to center is 8" - 12" (20-30cm). Kneepad depth range from back of heel is 3½" - 8" (9-20cm). Kneepad height range from footplate to top of pad is 17" - 22½" (43-57cm).



#PNG50044 Quad Grip Handle

Attaches to the actuator handle, for users with limited dexterity and function. Adjusts for various hand sizes.



#PNG50070 T-Style Handle

For people with limited range of motion. Makes the first few pumps to ascending easier and reaching to descend more accessible. Adds additional 6" (15cm) in each direction.



#P80829 Upper Body Support Strap

D-ring adjustment, provides additional upper posterior support while lifting and standing. Pad size is 5.5" x 24" (14x61cm). Adjustment range from 33" - 70" (84-178cm).



P82062 Grab Handles

Handles attach to the lifting arms to give the user a place to hold onto during the standing transition.



P82160 Transport Handle

Attaches to center post of unit to provide the caregiver with an easy way to move the unit from room to room.



P82188 Swing-Out Legs

With the simple shift of the lever, the caregiver can swing in/out the legs of the stander for wider wheelchair access. For people with wheelchair width greater than 22½" (58cm) and up to 38" (97cm) to access stander. Includes four locking casters.