



Retail Price List & Order Form

USA Effective May 1, 2021 Date: _____

Bantam Extra Small User Specifications

Approximate Height Range 28"-40"
 Weight Limit 50 lbs.
 Seat Depth Range (seat to back surface) 7"-12"
 Seat to Foot Plate Range (from seat pivot) 4.5"-16.5"

Bill To (EasyStand Supplier Only)

Purchase Order #: _____

Quote #: _____

Account #: _____

Ordered By: _____

Name of ATP (Required): _____

Supplier Company: _____

Address: _____

City, State, Zip: _____

Phone: _____

Fax: _____

Email: _____

Payer Source: _____

Ship To

Name: _____

Address: _____

Address: _____

City, State: _____

Zip: _____

Phone: _____

Tag For: _____

Email form to webquotes@easystand.com

Bantam Extra Small Packages

Select one Pattern and one Color

Hygienic Cover Patterns (select one)

- Rainbow Unicorn
- Sea Creatures
- Fancy Pants Elephants
- Space Race DEFAULT
- All Smiles Animals

Frame Color (select one)

- Green DEFAULT
- Purple
- White
- Orange

PK101 Minimum Support Package.....\$3,794

The minimum support package includes: (no substitutions)

Supine, Easy-Adjust Seat Depth, Black Molded Tray, Velcro® Positioning Belt, Head Support, Planar Seat & Back, Gas Spring Lift with Foot Pedal, Knee Pads, 5" Front Wheels, 3" Rear Locking Swivel Casters, Chest Strap, Foot Plates.

Foot Holders (select one, pair)

- 6"Lx3"W 7.75"Lx3.25"W DEFAULT 9.75"Lx4"W

Knee Pads (select one, inside width)

- 2.5" 3.25" DEFAULT 4.25"

Planar Back (select one, height from seat)

- 11"-13"H Range DEFAULT 15"-17"H Range

Velcro® Positioning Belt (select one, hip circumference)

- up to 23" DEFAULT up to 33"

Head Support (select one)

- 5"Hx8"W 6"Hx10"W 5"Hx10"W-Form to Fit DEFAULT

PK102 Moderate Support Package.....\$6,102

The moderate support package includes: (no substitutions)

Supine, Easy-Adjust Seat Depth, Black Molded Swing-Away Tray, Forearm Wings, Multi-Adjustable Foot Plates, Secure Foot Straps, Hygienic Covers for (Seat, Back, Head Support), Swing-Away Knee Pads, Form to Fit Seat, Form to Fit Back, Velcro® Positioning Belt, X-Style Chest Vest, High Mount Chest Vest Bracket, Head Support, Gas Spring Lift with Foot Pedal, Front Wheels/Casters, 3" Rear Locking Swivel Casters.

Front Wheel/Casters (select one)

- 5" Front Wheels DEFAULT Front Swivel Casters

Foot Holders (select one, pair)

- 6"Lx3"W 7.75"Lx3.25"W DEFAULT 9.75"Lx4"W

Secure Foot Straps (select one, length over top of foot, two pair)

- 10"L DEFAULT 15"L

Knee Pads (select one, inside width)

- 2.5" 3.25" DEFAULT 4.25"

From to Fit Back (select one, height from seat)

- 11"-13"H Range DEFAULT 15"-17"H Range

Velcro® Positioning Belt (select one, hip circumference)

- up to 23" DEFAULT up to 33"

X-Style Chest Vest (select one)

- 9.5"Lx9"W DEFAULT 11"Lx9.5"W
- Non-Stretch 9.5"Lx9"W Non-Stretch 11"Lx9.5"W

Head Support (select one)

- 5"Hx8"W 6"Hx10"W 5"Hx10"W-Form to Fit DEFAULT

PK103 Maximum Support Package.....\$6,740

The maximum support package includes: (no substitutions)

Supine, Easy-Adjust Seat Depth, Dual Control, Gas Spring Lockout, Swing-Away Tray, Forearm Wings, Multi-Adjustable Foot Plates, Secure Foot Straps, Hygienic Covers for (Seat, Back, Head Support), Swing-Away Knee Pads, Hip Supports, Calf Pad, Lateral Supports, Velcro® Positioning Belt, X-Style Chest Vest, High Mount Chest Vest Bracket, Head Support, Push Handle, Gas Spring Lift with Foot Pedal, Planar Seat & Back, Front Wheels/Casters, 3" Rear Locking Swivel Casters.

Front Wheel/Casters (select one)

- 5" Front Wheels DEFAULT Front Swivel Casters

Swing Away Tray (select one)

- Black Molded 11"Lx19"W Clear Tray-11"Lx19"W DEFAULT

Foot Holders (select one, pair)

- 6"Lx3"W 7.75"Lx3.25"W DEFAULT 9.75"Lx4"W

Secure Foot Straps (select one, length over top of foot, two pair)

- 10"L DEFAULT 15"L

Knee Pads (select one, inside width)

- 2.5" 3.25" DEFAULT 4.25"

Planar Back (select one, height from seat)

- 11"-13"H Range DEFAULT 15"-17"H Range

Lateral Supports (select one)

- 5"-9"W Range (select pad style below)
 - Flat Pads DEFAULT Curved Pads
- 8"-12"W Range (select pad style below)
 - Flat Pads DEFAULT Curved Pads

Velcro® Positioning Belt (select one, hip circumference)

- up to 23" DEFAULT up to 33"

X-Style Chest Vest (select one)

- 9.5"Lx9"W DEFAULT 11"Lx9.5"W
- Non-Stretch 9.5"Lx9"W Non-Stretch 11"Lx9.5"W

Head Support (select one)

- 5"Hx8"W 6"Hx10"W 5"Hx10"W-Form to Fit DEFAULT

PK104 Moderate Support Package (Mobile)....\$5,723

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Mobile, Mobile Tray, Multi-Adjustable Foot Plates, Secure Foot Straps, Hygienic Covers for (Seat, Back), Swing-Away Knee Pads, Hip Supports, Velcro® Positioning Belt, Chest Strap, X-Style Chest Vest, High Mount Chest Vest Bracket, Gas Spring Lift with Foot Pedal, Planar Seat & Back, 3" Rear Locking Casters, Chest Strap.

Foot Holders (select one, pair)

- 6"Lx3"W 7.75"Lx3.25"W DEFAULT 9.75"Lx4"W

Secure Foot Straps (select one, length over top of foot, two pair)

- 10"L DEFAULT 15"L

Knee Pads (select one, inside width)

- 2.5" 3.25" DEFAULT 4.25"

Planar Back (select one, height from seat)

- 11"-13"H Range DEFAULT 15"-17"H Range

Velcro® Positioning Belt (select one, hip circumference)

- up to 23" DEFAULT up to 33"

X-Style Chest Vest (select one)

- 9.5"Lx9"W DEFAULT 11"Lx9.5"W
- Non-Stretch 9.5"Lx9"W Non-Stretch 11"Lx9.5"W

AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

easystand.com



Retail Price List & Order Form

USA Effective May 1, 2021 Date: _____

Bantam Small User Specifications

Approximate Height Range 36"-54"
 Weight Limit 100 lbs.
 Seat Depth Range (seat to back surface) 11"-16"
 Seat to Foot Plate Range (from seat pivot) 4.5"-16.5"

Bill To (EasyStand Supplier Only)

Purchase Order #: _____

Quote #: _____

Account #: _____

Ordered By: _____

Name of ATP (Required): _____

Supplier Company: _____

Address: _____

City, State, Zip: _____

Phone: _____

Fax: _____

Email: _____

Payer Source: _____

Ship To

Name: _____

Address: _____

Address: _____

City, State: _____

Zip: _____

Phone: _____

Tag For: _____

Email form to webquotes@easystand.com

AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

easystand.com

Bantam Small Packages

Select one Pattern and one Color

Hygienic Cover Patterns (select one)

- Rainbow Unicorn
- Sea Creatures
- Fancy Pants Elephants
- Space Race DEFAULT
- All Smiles Animals

Frame Color (select one)

- Green DEFAULT
- Purple
- White
- Orange

PK105 Minimum Support Package.....\$3,794

The minimum support package includes: (no substitutions)

Supine, Easy-Adjust Seat Depth, Black Molded Tray, Velcro® Positioning Belt, Head Support, Planar Seat & Back, Gas Spring Lift with Foot Pedal, Knee Pads, 5" Front Wheels, 3" Rear Locking Swivel Casters, Chest Strap, Foot Plates.

Foot Holders (select one, pair)

- 7.75"Lx3.25"W 9.75"Lx4"W DEFAULT 11.75"Lx5"W

Knee Pads (select one, inside width)

- 3.25" 4.25" DEFAULT 5"

Planar Back (select one, height from seat)

- 14"-16"H Range DEFAULT 18"-20"H Range

Velcro® Positioning Belt (select one, hip circumference)

- up to 29" DEFAULT up to 39"

Head Support (select one)

- 5"Hx8"W 6"Hx10"W 5"Hx10"W-Form to Fit DEFAULT

PK106 Moderate Support Package.....\$6,102

The moderate support package includes: (no substitutions)

Supine, Easy-Adjust Seat Depth, Black Molded Swing-Away Tray, Forearm Wings, Multi-Adjustable Foot Plates, Secure Foot Straps, Hygienic Covers for (Seat, Back, Head Support), Swing-Away Knee Pads, Form to Fit Seat, Form to Fit Back, Velcro® Positioning Belt, X-Style Chest Vest, High Mount Chest Vest Bracket, Head Support, Gas Spring Lift with Foot Pedal, Front Wheels/Casters, 3" Rear Locking Swivel Casters.

Front Wheel/Casters (select one)

- 5" Front Wheels DEFAULT Front Swivel Casters

Foot Holders (select one, pair)

- 7.75"Lx3.25"W 9.75"Lx4"W DEFAULT 11.75"Lx5"W

Secure Foot Straps (select one, length over top of foot, two pair)

- 10"L DEFAULT 15"L

Knee Pads (select one, inside width)

- 3.25" 4.25" DEFAULT 5"

From to Fit Back (select one, height from seat)

- 12.5"-14.5"H Range DEFAULT 16.5"-18.5"H Range

Velcro® Positioning Belt (select one, hip circumference)

- up to 29" DEFAULT up to 39"

X-Style Chest Vest (select one)

- 9.5"Lx9"W DEFAULT 11"Lx9.5"W
- Non-Stretch 9.5"Lx9"W Non-Stretch 11"Lx9.5"W

Head Support (select one)

- 5"Hx8"W 6"Hx10"W 5"Hx10"W-Form to Fit DEFAULT

PK107 Maximum Support Package.....\$6,740

The maximum support package includes: (no substitutions)

Supine, Easy-Adjust Seat Depth, Dual Control, Gas Spring Lockout, Swing-Away Tray, Forearm Wings, Multi-Adjustable Foot Plates, Secure Foot Straps, Hygienic Covers for (Seat, Back, Head Support), Swing-Away Knee Pads, Hip Supports, Calf Pad, Lateral Supports, Velcro® Positioning Belt, X-Style Chest Vest, High Mount Chest Vest Bracket, Head Support, Push Handle, Gas Spring Lift with Foot Pedal, Planar Seat & Back, Front Wheels/Casters, 3" Rear Locking Swivel Casters.

Front Wheel/Casters (select one)

- 5" Front Wheels DEFAULT Front Swivel Casters

Swing Away Tray (select one)

- Black Molded 11"Lx19"W Clear Tray-11"Lx19"W DEFAULT

Knee Pads (select one, inside width)

- 3.25" 4.25" DEFAULT 5"

Foot Holders (select one, pair)

- 7.75"Lx3.25"W 9.75"Lx4"W DEFAULT 11.75"Lx5"W

Secure Foot Straps (select one, length over top of foot, two pair)

- 10"L DEFAULT 15"L

Planar Back (select one, height from seat)

- 14"-16"H Range DEFAULT 18"-20"H Range

Lateral Supports (select one)

- 5"-9"W Range (select pad style below)
 - Flat Pads DEFAULT Curved Pads
- 8"-12"W Range (select pad style below)
 - Flat Pads DEFAULT Curved Pads

Velcro® Positioning Belt (select one, hip circumference)

- up to 29" DEFAULT up to 39"

X-Style Chest Vest (select one)

- 9.5"Lx9"W DEFAULT 11"Lx9.5"W
- Non-Stretch 9.5"Lx9"W Non-Stretch 11"Lx9.5"W

Head Support (select one)

- 5"Hx8"W 6"Hx10"W 5"Hx10"W-Form to Fit DEFAULT

PK108 Moderate Support Package (Mobile)....\$5,723

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Mobile, Mobile Tray, Multi-Adjustable Foot Plates, Secure Foot Straps, Hygienic Covers for (Seat, Back), Swing-Away Knee Pads, Hip Supports, Velcro® Positioning Belt, Chest Strap, X-Style Chest Vest, High Mount Chest Vest Bracket, Gas Spring Lift with Foot Pedal, Planar Seat & Back, 3" Rear Locking Casters, Chest Strap.

Knee Pads (select one, inside width)

- 3.25" 4.25" DEFAULT 5"

Foot Holders (select one, pair)

- 7.75"Lx3.25"W 9.75"Lx4"W DEFAULT 11.75"Lx5"W

Secure Foot Straps (select one, length over top of foot, two pair)

- 10"L DEFAULT 15"L

Planar Back (select one, height from seat)

- 14"-16"H Range DEFAULT 18"-20"H Range

Velcro® Positioning Belt (select one, hip circumference)

- up to 29" DEFAULT up to 39"

X-Style Chest Vest (select one)

- 9.5"Lx9"W DEFAULT 11"Lx9.5"W
- Non-Stretch 9.5"Lx9"W Non-Stretch 11"Lx9.5"W



Retail Price List & Order Form

USA Effective May 1, 2021 Date: _____

Bantam Medium User Specifications

Approximate Height Range	4'0"-5'6"
Weight Limit	200 lbs.
Seat Depth Range (seat to back surface)	15"-20"
Seat to Foot Plate Range (from seat pivot)	7"-18.5"

Bill To (EasyStand Supplier Only)

Purchase Order #: _____

Quote #: _____

Account #: _____

Ordered By: _____

Name of ATP (Required): _____

Supplier Company: _____

Address: _____

City, State, Zip: _____

Phone: _____

Fax: _____

Email: _____

Payer Source: _____

Ship To

Name: _____

Address: _____

Address: _____

City, State: _____

Zip: _____

Phone: _____

Tax For: _____

Email form to webquotes@easystand.com

Bantam Medium Packages

PK109 Minimum Support Package.....\$5,723

The minimum support package includes: (no substitutions)
Supine, Easy-Adjust Seat Depth, Black Molded Tray, Planar Back, Velcro® Positioning Belt, 5" Front Wheels, Slide-In Knee Pads, 5" Rear Locking Casters, Planar Seat, Height-Only Adjustable Foot Plates.

Knee Pads (select one, inside width)

- 4" 5.5" DEFAULT

Foot Holders (select one, pair)

- 7.75"Lx3.25"W 9.75"Lx4"W 11.75"Lx5"W DEFAULT

Planar Back (select one, height from seat)

- 16"H 20"H DEFAULT

Head Support (select one)

- 5"Hx8"W 6"Hx10"W
- 7"Hx14"W 5"Hx10"W-Form to Fit DEFAULT

PK110 Moderate Support Package.....\$7,110

The moderate support package includes: (no substitutions)
Supine, Easy-Adjust Seat Depth, Black Molded Swing-Away Tray, Large Contoured Chest Pad, Multi-Adjustable Foot Plates, Secure Foot Straps, Swing-Away Knee Pads, Hygienic Covers (Seat, Back), Contoured Seat, Contoured Back, Velcro® Positioning Belt, Chest Strap or Chest Vest with High Mount Vest Bracket, Head Support, Front Frame, 5" Rear Locking Casters.

Front Frame (select one)

- 5" Front Wheels DEFAULT 5" Front Swivel Casters

Knee Pads (select one, inside width)

- 4" 5.5" DEFAULT

Foot Holders (select one, pair)

- 7.75"Lx3.25"W 9.75"Lx4"W 11.75"Lx5"W DEFAULT

Secure Foot Straps (select one, length over top of foot, two pair)

- 10"L DEFAULT 15"L

Contoured Seat (select one)

- Narrow 11"-13"W Wide 12"-14"W DEFAULT

Contoured Back (select one, height from seat)

- 16"H 20"H DEFAULT

Chest Support (select one)

- Chest Strap (circumference up to 42", velcro with D-ring)
- X-Style 11"Lx9.5"W X-Style 14.5"Lx10.5"W DEFAULT
- Y-Style 17.5"Lx11.5"W Y-Style 20.5"Lx12.5"W
- Non-Stretch X-Style 11"Lx9.5"W Non-Stretch X-Style 14.5"Lx10.5"W
- Non-Stretch Y-Style 17.5"Lx11.5"W Non-Stretch Y-Style 20.5"Lx12.5"W

Head Support (select one)

- 5"Hx8"W 6"Hx10"W
- 7"Hx14"W 5"Hx10"W-Form to Fit DEFAULT

PK111 Maximum Support Package.....\$7,818

The maximum support package includes: (no substitutions)
Supine, Easy-Adjust Seat Depth, Swing-Away Tray, Large Contoured Chest Pad, Multi-Adjustable Foot Plates, Secure Foot Straps, Swing-Away Knee Pads, Hygienic Covers (Seat, Back), Hip Supports, Planar Back, Velcro® Positioning Belt, Chest Strap or Chest Vest with High Mount Vest Bracket, Lateral Supports, Head Support, Front Frame, 5" Rear Locking Casters, Planar Seat.

Front Frame (select one)

- 5" Front Wheels DEFAULT 5" Front Swivel Casters

Swing Away Tray (select one)

- Black Molded 21"Lx24"W Clear Tray-21"Lx24"W DEFAULT

Knee Pads (select one, inside width)

- 4" 5.5" DEFAULT

Foot Holders (select one, pair)

- 7.75"Lx3.25"W 9.75"Lx4"W 11.75"Lx5"W DEFAULT

Secure Foot Straps (select one, length over top of foot, two pair)

- 10"L DEFAULT 15"L

Planar Back (select one, height from seat)

- 16"H 20"H DEFAULT

Chest Support (select one)

- Chest Strap (circumference up to 42", velcro with D-ring)
- X-Style 11"Lx9.5"W X-Style 14.5"Lx10.5"W DEFAULT
- Y-Style 17.5"Lx11.5"W Y-Style 20.5"Lx12.5"W
- Non-Stretch X-Style 11"Lx9.5"W Non-Stretch X-Style 14.5"Lx10.5"W
- Non-Stretch Y-Style 17.5"Lx11.5"W Non-Stretch Y-Style 20.5"Lx12.5"W

Lateral Supports (select one)

- 9"-15.5"W Range (select pad style below)
 - Flat Pads Curved Pads
- With Elbow Stop 9"-14.5" Range DEFAULT

Head Support (select one)

- 5"Hx8"W 6"Hx10"W
- 7"H x 14"W 5"Hx10"W-Form to Fit DEFAULT

AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

easystand.com



Retail Price List & Order Form

USA Effective May 1, 2021 Date: _____

Evolv Medium User Specifications

Approximate Height Range	4'0"-5'6"
Weight Limit	200 lbs.
Seat Depth Range (seat to back surface)	14"-19"
Seat to Foot Plate Range (from seat pivot)	11.5"-18"

Bill To (EasyStand Supplier Only)

Purchase Order #: _____

Quote #: _____

Account #: _____

Ordered By: _____

Name of ATP (Required): _____

Supplier Company: _____

Address: _____

City, State, Zip: _____

Phone: _____

Fax: _____

Email: _____

Payer Source: _____

Ship To

Name: _____

Address: _____

Address: _____

City, State: _____

Zip: _____

Phone: _____

Tag For: _____

Email form to webquotes@easystand.com

AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

easystand.com

Evolv Medium Packages

PK201 Minimum Support Package.....\$4,307

The minimum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Removable Back, Velcro® Positioning Belt, Chest Strap, Black Molded Tray with Chest Pad, Manual Hydraulic Actuator with Handle, Flip-Up Knee Pad, Planar Seat, Multi-Adjustable Foot Plates, 5" Front Wheels, 5" Rear Locking Casters.

Actuator Handle Location (select one)

Right DEFAULT Left

Removable Back (select one)

Flat Back-15"H Flat Back-19"H DEFAULT
 Contoured Back-15"H Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

26"-55" DEFAULT 38"-67"

Chest Strap (select one, chest measurement)

30"-46" DEFAULT 34"-54" 46"-66"

PK202 Moderate Support Package.....\$6,322

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Swing-Away Front, Adjustable/Removable Actuator Handle, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads for Swing-Away, Removable Back, Velcro® Positioning Belt, Chest Strap, Accessories Mounting Bracket, Lateral Supports, Planar Seat, Multi-Adjustable Foot Plates, Black Molded Tray, 5" Front Wheels, 5" Rear Locking Casters.

Actuator Handle Location (select one)

Right DEFAULT Left

Chest Pad (select one)

Non-Adjustable DEFAULT Adjustable

Secure Foot Straps (select one, length over top of foot, two pair)

15"L DEFAULT 21"L

Hip Supports (select one range)

9"-15" DEFAULT 11"-17" 16"-22"

Removable Back (select one)

Flat Back-15"H Flat Back-19"H DEFAULT
 Contoured Back-15"H Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

26"-55" DEFAULT 38"-67"

Chest Strap (select one, chest measurement)

30"-46" DEFAULT 34"-54" 46"-66"

Lateral Supports (select one)

8"-16.5"W Range (select pad style below)
 Flat Pads DEFAULT Curved Pads
 11"-19.5"W Range (select pad style below)
 Flat Pads Curved Pads

PK203 Maximum Support Package.....\$7,989

The maximum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Swing-Away Front for Shadow Tray, Adjustable/Removable Actuator Handle, Shadow Tray, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads for Swing-Away, Removable Back, Positioning Belt, Accessories Mounting Bracket, Chest Strap or Chest Vest with High Mount Vest Bracket, Lateral Supports, Head Support, Tray, Multi-Adjustable Foot Plates, Planar Seat, Front Wheels/Casters, 5" Rear Locking Casters.

Front Wheel/Casters (select one)

5" Front Wheels DEFAULT Front Swivel Casters

Actuator Handle Location (select one)

Right DEFAULT Left

Shadow Tray (select one)

Black Molded Adjustable Shadow Tray (7"-12"H from seat, 4"-20"D from back)
 Clear Adjustable Shadow Tray (6"-11"H from seat, 4"-20"D from back) DEFAULT

Chest Pad (select one)

Non-Adjustable DEFAULT Adjustable

Secure Foot Straps (select one, length over top of foot, two pair)

15"L DEFAULT 21"L

Hip Supports (select one range)

9"-15" DEFAULT 11"-17" 16"-22"

Removable Back (select one)

Flat Back-15"H Flat Back-19"H
 Contoured Back-15"H Contoured Back-19"H DEFAULT

Velcro® Positioning Belt (select one, hip circumference)

26"-55" DEFAULT 38"-67"

Chest Support (select one)

Chest Strap select one chest measurement: 30"-46" 34"-54" 46"-66"
 X-Style 11"Lx9.5"W X-Style 14.5"Lx10.5"W DEFAULT
 Y-Style 17.5"Lx11.5"W Y-Style 20.5"Lx12.5"W
 Non-Stretch X-Style 11"Lx9.5"W Non-Stretch X-Style 14.5"Lx10.5"W
 Non-Stretch Y-Style 17.5"Lx11.5"W Non-Stretch Y-Style 20.5"Lx12.5"W

Lateral Supports (select one)

8"-16.5"W Range (select pad style below)
 Flat Pads DEFAULT Curved Pads
 11"-19.5"W Range (select pad style below)
 Flat Pads Curved Pads

Head Support (select one)

Short DEFAULT Tall

PK204 Minimum Support Package (Mobile)....\$6,403

The minimum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Mobile, Clear Angle Adjustable Tray for Mobile, Large Contoured Chest Pad, Removable Back, Velcro® Positioning Belt, Chest Strap, Manual Hydraulic Actuator with Handle, Planar Seat, Flip-Up Knee Pads, Multi-Adjustable Foot Plates, 5" Rear Locking Casters.

Actuator Handle Location (select one)

Right DEFAULT Left

Chest Pad (select one)

Non-Adjustable DEFAULT Adjustable

Removable Back (select one)

Flat Back-15"H Flat Back-19"H DEFAULT
 Contoured Back-15"H Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

26"-55" DEFAULT 38"-67"

Chest Strap (select one, chest measurement)

30"-46" DEFAULT 34"-54" 46"-66"



Retail Price List & Order Form

USA Effective May 1, 2021 Date: _____

Evolv Large User Specifications

Approximate Height Range	5'0"-6'2"
Weight Limit	280 lbs.
Seat Depth Range (seat to back surface)	18"-23"
Seat to Foot Plate Range (from seat pivot)	14"-21"

Bill To (EasyStand Supplier Only)

Purchase Order #: _____

Quote #: _____

Account #: _____

Ordered By: _____

Name of ATP (Required): _____

Supplier Company: _____

Address: _____

City, State, Zip: _____

Phone: _____

Fax: _____

Email: _____

Payer Source: _____

Ship To

Name: _____

Address: _____

Address: _____

City, State: _____

Zip: _____

Phone: _____

Tax For: _____

Email form to webquotes@easystand.com

Evolv Large Packages

PK205 Minimum Support Package.....\$4,307

The minimum support package includes: (no substitutions)
Easy-Adjust Seat Depth, Removable Back, Velcro® Positioning Belt, Chest Strap, Black Molded Tray with Chest Pad, Manual Hydraulic Actuator with Handle, Flip-Up Knee Pad, Planar Seat, Multi-Adjustable Foot Plates, 5" Front Wheels, 5" Rear Locking Casters.

Actuator Handle Location (select one)

- Right DEFAULT Left

Removable Back (select one)

- Flat Back-15"H Flat Back-19"H DEFAULT
 Contoured Back-15"H Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

- 26"-55" DEFAULT 38"-67"

Chest Strap (select one, chest measurement)

- 30"-46" 34"-54" DEFAULT 46"-66"

PK206 Moderate Support Package.....\$6,322

The moderate support package includes: (no substitutions)
Easy-Adjust Seat Depth, Swing-Away Front, Adjustable/Removable Actuator Handle, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads for Swing-Away, Removable Back, Velcro® Positioning Belt, Chest Strap, Accessories Mounting Bracket, Lateral Supports, Planar Seat, Multi-Adjustable Foot Plates, Black Molded Tray, 5" Front Wheels, 5" Rear Locking Casters.

Actuator Handle Location (select one)

- Right DEFAULT Left

Chest Pad (select one)

- Non-Adjustable DEFAULT Adjustable

Secure Foot Straps (select one, length over top of foot, two pair)

- 15"L DEFAULT 21"L

Hip Supports (select one range)

- 10"-15" 12"-17" DEFAULT 17"-23"

Removable Back (select one)

- Flat Back-15"H Flat Back-19"H DEFAULT
 Contoured Back-15"H Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

- 26"-55" DEFAULT 38"-67"

Chest Strap (select one, chest measurement)

- 30"-46" 34"-54" DEFAULT 46"-66"

Lateral Supports (select one)

- 8"-16.5"W Range (select pad style below)
 Flat Pads Curved Pads
 11"-19.5"W Range (select pad style below)
 Flat Pads DEFAULT Curved Pads

PK207 Maximum Support Package.....\$7,989

The maximum support package includes: (no substitutions)
Easy-Adjust Seat Depth, Swing-Away Front for Shadow Tray, Adjustable/Removable Actuator Handle, Shadow Tray, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads for Swing-Away, Removable Back, Positioning Belt, Accessories Mounting Bracket, Chest Strap or Chest Vest with High Mount Vest Bracket, Lateral Supports, Head Support, Tray, Multi-Adjustable Foot Plates, Planar Seat, Front Wheels/Casters, 5" Rear Locking Casters.

Front Wheel/Casters (select one)

- 5" Front Wheels DEFAULT Front Swivel Casters

Actuator Handle Location (select one)

- Right DEFAULT Left

Shadow Tray (select one)

- Black Molded Adjustable Shadow Tray (7"-12"H from seat, 4"-20"D from back)
 Clear Adjustable Shadow Tray (6"-11"H from seat, 4"-20"D from back) DEFAULT

Chest Pad (select one)

- Non-Adjustable DEFAULT Adjustable

Secure Foot Straps (select one, length over top of foot, two pair)

- 15"L DEFAULT 21"L

Hip Supports (select one range)

- 10"-15" 12"-17" DEFAULT 17"-23"

Removable Back (select one)

- Flat Back-15"H Flat Back-19"H
 Contoured Back-15"H Contoured Back-19"H DEFAULT

Velcro® Positioning Belt (select one, hip circumference)

- 26"-55" DEFAULT 38"-67"

Chest Support (select one)

- Chest Strap select one chest measurement: 30"-46" 34"-54" 46"-66"
 X-Style 11"Lx9.5"W X-Style 14.5"Lx10.5"W DEFAULT
 Y-Style 17.5"Lx11.5"W Y-Style 20.5"Lx12.5"W
 Non-Stretch X-Style 11"Lx9.5"W Non-Stretch X-Style 14.5"Lx10.5"W
 Non-Stretch Y-Style 17.5"Lx11.5"W Non-Stretch Y-Style 20.5"Lx12.5"W

Lateral Supports (select one)

- 8"-16.5"W Range (select pad style below)
 Flat Pads Curved Pads
 11"-19.5"W Range (select pad style below)
 Flat Pads DEFAULT Curved Pads

Head Support (select one)

- Short Tall DEFAULT

PK208 Minimum Support Package (Mobile)....\$6,403

The minimum support package includes: (no substitutions)
Easy-Adjust Seat Depth, Mobile, Clear Angle Adjustable Tray for Mobile, Large Contoured Chest Pad, Removable Back, Velcro® Positioning Belt, Chest Strap, Manual Hydraulic Actuator with Handle, Planar Seat, Flip-Up Knee Pads, Multi-Adjustable Foot Plates, 5" Rear Locking Casters.

Actuator Handle Location (select one)

- Right DEFAULT Left

Chest Pad (select one)

- Non-Adjustable DEFAULT Adjustable

Removable Back (select one)

- Flat Back-15"H Flat Back-19"H DEFAULT
 Contoured Back-15"H Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

- 26"-55" DEFAULT 38"-67"

Chest Strap (select one, chest measurement)

- 30"-46" 34"-54" DEFAULT 46"-66"

AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

easystand.com



Evolv XT Packages

Retail Price List & Order Form

USA Effective May 1, 2021 Date: _____

Evolv XT User Specifications

Approximate Height Range	6'0"-6'10"
Weight Limit	350 lbs.
Seat Depth Range (seat to back surface)	20"-25"
Seat to Foot Plate Range (from seat pivot)	17"-23"

Bill To (EasyStand Supplier Only)

Purchase Order #: _____

Quote #: _____

Account #: _____

Ordered By: _____

Name of ATP (Required): _____

Supplier Company: _____

Address: _____

City, State, Zip: _____

Phone: _____

Fax: _____

Email: _____

Payer Source: _____

Ship To

Name: _____

Address: _____

Address: _____

City, State: _____

Zip: _____

Phone: _____

Tax For: _____

Email form to webquotes@easystand.com

AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

easystand.com

PK209 Minimum Support Package.....\$5,579

The minimum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Removable Back, Velcro® Positioning Belt, Chest Strap, Black Molded Tray with Chest Pad, Manual Hydraulic Actuator with Handle, Flip-Up Knee Pad, Planar Seat, Multi-Adjustable Foot Plates, 5" Front Wheels, 5" Rear Locking Casters.

Actuator Handle Location (select one)

Right DEFAULT Left

Removable Back (select one)

Flat Back-15"H Flat Back-19"H DEFAULT
 Contoured Back-15"H Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

26"-55" 38"-67" DEFAULT

Chest Strap (select one, chest measurement)

30"-46" 34"-54" 46"-66" DEFAULT

PK210 Moderate Support Package.....\$7,128

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Adjustable/Removable Actuator Handle, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads, Removable Back, Velcro® Positioning Belt, Chest Strap, Accessories Mounting Bracket, Lateral Supports, Planar Seat, Multi-Adjustable Foot Plates, Black Molded Tray, 5" Front Wheels, 5" Rear Locking Casters.

Actuator Handle Location (select one)

Right DEFAULT Left

Chest Pad (select one)

Non-Adjustable DEFAULT Adjustable

Secure Foot Straps (select one, length over top of foot, two pair)

15"L DEFAULT 21"L

Hip Supports (select one range)

10"-15" 12"-17" 17"-23" DEFAULT

Removable Back (select one)

Flat Back-15"H Flat Back-19"H DEFAULT
 Contoured Back-15"H Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

26"-55" 38"-67" DEFAULT

Chest Strap (select one, chest measurement)

30"-46" 34"-54" 46"-66" DEFAULT

Lateral Supports (select one)

8"-16.5"W Range (select pad style below)
 Flat Pads Curved Pads
 11"-19.5"W Range (select pad style below)
 Flat Pads DEFAULT Curved Pads

PK211 Maximum Support Package.....\$8,851

The maximum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Adjustable/Removable Actuator Handle, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads, Removable Back, Velcro® Positioning Belt, Accessories Mounting Bracket, Chest Strap or Chest Vest with High Mount Vest Bracket, Lateral Supports, Head Support, Tray, Multi-Adjustable Foot Plates, Planar Seat, 5" Front Wheels/ Casters, 5" Rear Locking Casters.

Front Wheel/Casters (select one)

5" Front Wheels DEFAULT Front Swivel Casters

Actuator Handle Location (select one)

Right DEFAULT Left

Shadow Tray (select one)

Black Molded Adjustable Shadow Tray (7"-12"H from seat, 4"-20"D from back)
 Clear Adjustable Shadow Tray (6"-11"H from seat, 4"-20"D from back) DEFAULT

Chest Pad (select one)

Non-Adjustable DEFAULT Adjustable

Secure Foot Straps (select one, length over top of foot, two pair)

15"L DEFAULT 21"L

Hip Supports (select one range)

10"-15" 17"-17" 17"-23" DEFAULT

Removable Back (select one)

Flat Back-15"H Flat Back-19"H
 Contoured Back-15"H Contoured Back-19"H DEFAULT

Velcro® Positioning Belt (select one, hip circumference)

26"-55" 38"-67" DEFAULT

Chest Support (select one)

Chest Strap select one chest measurement: 30"-46" 34"-54" 46"-66"
 X-Style 11"Lx9.5"W X-Style 14.5"Lx10.5"W DEFAULT
 Y-Style 17.5"Lx11.5"W Y-Style 20.5"Lx12.5"W
 Non-Stretch X-Style 11"Lx9.5"W Non-Stretch X-Style 14.5"Lx10.5"W
 Non-Stretch Y-Style 17.5"Lx11.5"W Non-Stretch Y-Style 20.5"Lx12.5"W

Lateral Supports (select one)

8"-16.5"W Range (select pad style below)
 Flat Pads Curved Pads
 11"-19.5"W Range (select pad style below)
 Flat Pads DEFAULT Curved Pads

Head Support (select one)

Short Tall DEFAULT

PK212 Minimum Support Package (Mobile)....\$7,408

The minimum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Mobile, Clear Angle Adjustable Tray for Mobile, Large Contoured Chest Pad, Removable Back, Velcro® Positioning Belt, Chest Strap, Manual Hydraulic Actuator with Handle, Planar Seat, Flip-Up Knee Pads, Multi-Adjustable Foot Plates, 5" Rear Locking Casters.

Actuator Handle Location (select one)

Right DEFAULT Left

Chest Pad (select one)

Non-Adjustable DEFAULT Adjustable

Removable Back (select one)

Flat Back-15"H Flat Back-19"H DEFAULT
 Contoured Back-15"H Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

26"-55" 38"-67" DEFAULT

Chest Strap (select one, chest measurement)

30"-46" 34"-54" DEFAULT 46"-66"



Retail Price List & Order Form

USA Effective May 1, 2021 Date: _____

Glider User Specifications

This price list is for both the Glider Medium and the Glider Large. Specifications for each are listed above the pricing.

Bill To (EasyStand Supplier Only)

Purchase Order #: _____

Quote #: _____

Account #: _____

Ordered By: _____

Name of ATP (Required): _____

Supplier Company: _____

Address: _____

City, State, Zip: _____

Phone: _____

Fax: _____

Email: _____

Payer Source: _____

Ship To

Name: _____

Address: _____

Address: _____

City, State: _____

Zip: _____

Phone: _____

Tag For: _____

Email form to webquotes@easystand.com

AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

easystand.com

Glider Medium & Large Packages

Glider Medium User Specifications

Approximate Height Range	4'0"-5'6"
Weight Limit	200 lbs.
Seat Depth Range (center of knee pad to flat back)	20"-25"
Seat to Foot Plate Range (from seat pivot)	13.5"-16.5"

Glider Large User Specifications

Approximate Height Range	5'0"-6'2"
Weight Limit	280 lbs.
Seat Depth Range (center of knee pad to flat back)	24"-29"
Seat to Foot Plate Range (from seat pivot)	15.5"-18.5"

EasyStand Glider Medium

PK213 Moderate Support Package.....\$8,911

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Adjustable/Removable Actuator Handle, Secure Foot Straps, Hip Supports, Glide Handle Extensions, Removable Back, Velcro® Positioning Belt, Chest Strap, Full Range Seat, Positioning Bar, 5" Front Wheels, 5" Rear Locking Casters, Tray & Chest Pad, Flip-Up Knee Pads, Adjustable Resistance Cylinders, Adjustable Foot Plates.

Actuator Handle Location (select one)

- Right DEFAULT Left

Secure Foot Straps (select one, length over top of foot, two pair)

- 15"L DEFAULT 21"L

Hip Supports (select one range)

- 9"-15" DEFAULT 11"-17"

Removable Back (select one)

- Flat Back-15"H Flat Back-19"H DEFAULT
 Contoured Back-15"H Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

- 26"-55" DEFAULT 38"-67"

Chest Strap (select one, chest measurement)

- 30"-46" 34"-54" DEFAULT 46"-66"

EasyStand Glider Large

PK214 Moderate Support Package.....\$8,764

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Adjustable/Removable Actuator Handle, Secure Foot Straps, Hip Supports, Removable Back, Velcro® Positioning Belt, Chest Strap, Full Range Seat, Positioning Bar, 5" Front Wheels, 5" Rear Locking Casters, Tray & Chest Pad, Flip-Up Knee Pads, Adjustable Resistance Cylinders, Glide Handles, Adjustable Foot Plates.

Actuator Handle Location (select one)

- Right DEFAULT Left

Secure Foot Straps (select one, length over top of foot, two pair)

- 15"L DEFAULT 21"L

Hip Supports (select one range)

- 10"-15" 12"-17" DEFAULT

Removable Back (select one)

- Flat Back-15"H Flat Back-19"H DEFAULT
 Contoured Back-15"H Contoured Back-19"H

Velcro® Positioning Belt (select one, hip circumference)

- 26"-55" 38"-67" DEFAULT

Chest Strap (select one, chest measurement)

- 30"-46" 34"-54" DEFAULT 46"-66"



Retail Price List & Order Form

USA Effective May 1, 2021 Date: _____

StrapStand User Specifications

Approximate Height Range 5'0"-6'5"
Weight Limit 350 lbs.

Bill To (EasyStand Supplier Only)

Purchase Order #: _____

Quote #: _____

Account #: _____

Ordered By: _____

Name of ATP (Required): _____

Supplier Company: _____

Address: _____

City, State, Zip: _____

Phone: _____

Fax: _____

Email: _____

Payer Source: _____

Ship To

Name: _____

Address: _____

Address: _____

City, State: _____

Zip: _____

Phone: _____

Tag For: _____

Email form to webquotes@easystand.com

AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

easystand.com

StrapStand Packages

PK301 Minimum Support Package.....\$3,986

The minimum support package includes: (no substitutions)

Swing-Out Legs, Adjustable Lifting Strap, Black Molded Tray, Four Locking Swivel Casters, Foot Plates, Knee Pad, Hydraulic Actuator with Handle.

Lifting Strap (select one)

- Adjustable Lifting Strap 10"Wx26"L
- Adjustable Lifting Strap 12"Wx29"L DEFAULT
- Adjustable Lifting Strap 14"Wx32"L
- Adjustable Sling Strap 9"Wx25"L
- Adjustable Sling Strap 12"Wx28"L
- Adjustable Sling Strap 14"Wx32"L

PK302 Moderate Support Package.....\$5,532

The moderate support package includes: (no substitutions)

Swing-Out Legs, Adjustable Lifting Strap, Black Molded Tray, Four Locking Swivel Casters. Multi-Adjustable Footplates, Foot Straps, Independent Knee Pads, Removable Actuator Handle, Upper Body Support Strap, Grab Handles, Transport Handle.

Lifting Strap (select one)

- Adjustable Lifting Strap 10"Wx26"L
- Adjustable Lifting Strap 12"Wx29"L DEFAULT
- Adjustable Lifting Strap 14"Wx32"L
- Adjustable Sling Strap 9"Wx25"L
- Adjustable Sling Strap 12"Wx28"L
- Adjustable Sling Strap 14"Wx32"L

Secure Foot Straps (select one, length over top of foot, two pair)

- 15"L
- 21"L DEFAULT

PK303 Maximum Support Package.....\$5,898

The maximum support package includes: (no substitutions)

Swing-Out Legs, Adjustable Lifting Strap, Black Molded Tray, Large Contoured Chest Pad, Padded Tray Cover, Multi-Adjustable Footplates, Secure Foot Straps, Independent Roho Knee Pads, Removable Actuator Handle, Upper Body Support Strap, Grab Handles, Transport Handle, Four Locking Swivel Casters.

Lifting Strap (select one)

- Adjustable Lifting Strap 10"Wx26"L
- Adjustable Lifting Strap 12"Wx29"L DEFAULT
- Adjustable Lifting Strap 14"Wx32"L
- Adjustable Sling Strap 9"Wx25"L
- Adjustable Sling Strap 12"Wx28"L
- Adjustable Sling Strap 14"Wx32"L

Secure Foot Straps (select one, length over top of foot, two pair)

- 15"L
- 21"L DEFAULT