



# Retail Order Form

USA Effective January 1, 2024 Date: \_\_\_\_\_

## Evolv XT User Specifications

|  |            |
|--|------------|
| Approximate Height Range                   | 6'0"-6'10" |
| Weight Limit                               | 350 lbs.   |
| Seat Depth Range (seat to back surface)    | 20"-25"    |
| Seat to Foot Plate Range (from seat pivot) | 17"-23"    |

### Bill To (EasyStand Supplier Only)

Purchase Order #: \_\_\_\_\_

Quote #: \_\_\_\_\_

Account #: \_\_\_\_\_

Ordered By: \_\_\_\_\_

Name of ATP (Required): \_\_\_\_\_

Supplier Company: \_\_\_\_\_

Address: \_\_\_\_\_

City, State, Zip: \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

Payer Source: \_\_\_\_\_

### Ship To

Name: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

City, State: \_\_\_\_\_

Zip: \_\_\_\_\_

Phone: \_\_\_\_\_

Tag For: \_\_\_\_\_

Email form to [webquotes@easystand.com](mailto:webquotes@easystand.com)

# AltimateMedical

Ph: 800.342.8968 or 507.697.6393  
Fax: 877.342.8968 or 507.697.6900  
[easystand.com](http://easystand.com)

# Evolv XT Packages

## PK209 Minimum Support Package.....\$7,218

The minimum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Removable Back, Velcro® Positioning Belt, Chest Strap, Black Molded Tray with Chest Pad, Manual Hydraulic Actuator with Handle, Flip-Up Knee Pad, Planar Seat, Multi-Adjustable Foot Plates, 5" Front Wheels, 5" Rear Locking Casters.

### Actuator Handle Location (select one)

- Right DEFAULT  Left

### Removable Back (select one)

- Flat Back-15"H  Flat Back-19"H DEFAULT  
 Contoured Back-15"H  Contoured Back-19"H

### Velcro® Positioning Belt (select one, hip circumference)

- 26"-55"  38"-67" DEFAULT

### Chest Strap (select one, chest measurement)

- 30"-46"  34"-54"  46"-66" DEFAULT

## PK210 Moderate Support Package.....\$9,239

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Adjustable/Removable Actuator Handle, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads, Removable Back, Velcro® Positioning Belt, Chest Strap, Accessories Mounting Bracket, Lateral Supports, Planar Seat, Multi-Adjustable Foot Plates, Black Molded Tray, 5" Front Wheels, 5" Rear Locking Casters.

### Actuator Handle Location (select one)

- Right DEFAULT  Left

### Chest Pad (select one)

- Non-Adjustable DEFAULT  Adjustable

### Secure Foot Straps (select one, length over top of foot, two pair)

- 15"L DEFAULT  21"L

### Hip Supports (select one range)

- 10"-15"  12"-17"  17"-23" DEFAULT

### Removable Back (select one)

- Flat Back-15"H  Flat Back-19"H DEFAULT  
 Contoured Back-15"H  Contoured Back-19"H

### Velcro® Positioning Belt (select one, hip circumference)

- 26"-55"  38"-67" DEFAULT

### Chest Strap (select one, chest measurement)

- 30"-46"  34"-54"  46"-66" DEFAULT

### Lateral Supports (select one)

- 8"-16.5"W Range (select pad style below)  
 Flat Pads  Curved Pads  
 11"-19.5"W Range (select pad style below)  
 Flat Pads DEFAULT  Curved Pads

## PK211 Maximum Support Package.....\$10,708

The maximum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Adjustable/Removable Actuator Handle, Large Contoured Chest Pad, Secure Foot Straps, Hip Supports, Independent Knee Pads, Removable Back, Velcro® Positioning Belt, Accessories Mounting Bracket, Chest Strap or Chest Vest with High Mount Vest Bracket, Lateral Supports, Head Support, Shadow Tray, Multi-Adjustable Foot Plates, Planar Seat, 5" Front Wheels/Casters, 5" Rear Locking Casters.

### Front Wheel/Casters (select one)

- 5" Front Wheels DEFAULT  Front Swivel Casters

### Actuator Handle Location (select one)

- Right DEFAULT  Left

### Shadow Tray (select one)

- Black Molded Adjustable Shadow Tray (7"-12"H from seat, 4"-20"D from back)  
 Clear Adjustable Shadow Tray (6"-11"H from seat, 4"-20"D from back) DEFAULT

### Chest Pad (select one)

- Non-Adjustable DEFAULT  Adjustable

### Secure Foot Straps (select one, length over top of foot, two pair)

- 15"L DEFAULT  21"L

### Hip Supports (select one range)

- 10"-15"  17"-17"  17"-23" DEFAULT

### Removable Back (select one)

- Flat Back-15"H  Flat Back-19"H  
 Contoured Back-15"H  Contoured Back-19"H DEFAULT

### Velcro® Positioning Belt (select one, hip circumference)

- 26"-55"  38"-67" DEFAULT

### Chest Support (select one)

- Chest Strap select one chest measurement:  30"-46"  34"-54"  46"-66"  
 X-Style 11"Lx9.5"W  X-Style 14.5"Lx10.5"W DEFAULT  
 Y-Style 17.5"Lx11.5"W  Y-Style 20.5"Lx12.5"W  
 Non-Stretch X-Style 11"Lx9.5"W  Non-Stretch X-Style 14.5"Lx10.5"W  
 Non-Stretch Y-Style 17.5"Lx11.5"W  Non-Stretch Y-Style 20.5"Lx12.5"W

### Lateral Supports (select one)

- 8"-16.5"W Range (select pad style below)  
 Flat Pads  Curved Pads  
 11"-19.5"W Range (select pad style below)  
 Flat Pads DEFAULT  Curved Pads

### Head Support (select one)

- Short  Tall DEFAULT

## PK212 Minimum Support Package (Mobile)....\$9,988

The minimum support package includes: (no substitutions)

Easy-Adjust Seat Depth, Mobile, Clear Angle Adjustable Tray for Mobile, Large Contoured Chest Pad, Removable Back, Velcro® Positioning Belt, Chest Strap, Manual Hydraulic Actuator with Handle, Planar Seat, Flip-Up Knee Pads, Multi-Adjustable Foot Plates, 5" Rear Locking Casters.

### Actuator Handle Location (select one)

- Right DEFAULT  Left

### Chest Pad (select one)

- Non-Adjustable DEFAULT  Adjustable

### Removable Back (select one)

- Flat Back-15"H  Flat Back-19"H DEFAULT  
 Contoured Back-15"H  Contoured Back-19"H

### Velcro® Positioning Belt (select one, hip circumference)

- 26"-55"  38"-67" DEFAULT

### Chest Strap (select one, chest measurement)

- 30"-46"  34"-54" DEFAULT  46"-66"