

Bill TO (Ultimate Medical Supplier Only) 

Purchase Order #:	_____
Quote #:	_____
Account #:	_____
Ordered By:	_____
Supplier Company:	_____
Address:	_____
City, State, Zip:	_____
Phone:	_____ Email: _____
Payer Source:	_____
Name of ATP (Required):	_____

Ship To 

Name:	_____
Address:	_____
City, State, Zip:	_____
Phone:	_____
Tag For:	_____

Glider Large User Specifications

Approximate Height Range	5'0"-6'2"
Weight Limit	280 lbs.
Seat Depth Range (center of knee pad to flat back)	24"-29"
Seat to Foot Plate Range (from seat pivot)	15.5"-18.5"

Base *Note: The DEFAULT items will be added if no selection is made.*

EasyStand® Glider Large	\$10,087.51
Price Includes:	
PNG50378 Glider Base	\$4,612.01
PNG50024 Glider Configuration	\$5,475.50

Lift Mechanisms-must select one

<input type="checkbox"/> PNG50009 Manual Hydraulic Actuator DEFAULT	N/C
Manual Hydraulic Actuator Handle Location	
<input type="checkbox"/> PNG50059 Right Side DEFAULT	N/C
<input type="checkbox"/> PNG50126 Left Side	N/C
Manual Hydraulic Actuator Handle Accessories	
<input type="checkbox"/> PNG50328 Adjustable/Removable Actuator Handle	\$317.42
<input type="checkbox"/> PNG50044 Quad Grip Handle Extension	\$213.86
<input type="checkbox"/> PNG50070 T-Style Handle Extension	\$213.86
<input type="checkbox"/> PNG50010 Pow'r Up Lift (battery powered)	\$1,701.43

Tray-must select one

<input type="checkbox"/> PNG50370 Glider Tray DEFAULT	\$465.37
--	----------

Foot Plate Accessories

<input type="checkbox"/> PNG30030 Foot Straps-15"L (length over top of foot, pair)	\$67.25
<input type="checkbox"/> PNG50505 Foot Straps-21"L (length over top of foot, pair)	\$67.25
<input type="checkbox"/> PNG30031 Secure Foot Straps-15"L (length over top of foot, two pair)	\$169.47
<input type="checkbox"/> PNG50506 Secure Foot Straps-21"L (length over top of foot, two pair)	\$169.47

Seats-must select one

<input type="checkbox"/> PNG50352-1 Full Range Seat DEFAULT	\$815.07
<input type="checkbox"/> PNG50353 Full Range Seat with Roho® Insert	\$1,058.52

Seat Accessories

<input type="checkbox"/> PNG50499 Easy-Adjust Seat Depth- CANNOT be ordered off-unit (roller bearing allows for smoother seat depth adjustment)	\$150.64
--	----------

Pelvic Support

<input type="checkbox"/> PNG50345 Hip Supports (10"-15" range, pair)	\$505.72
<input type="checkbox"/> PNG50374 Hip Supports (12"-17" range, pair)	\$505.72
<input type="checkbox"/> PNG30029 Velcro® Positioning Belt-26"-55" (back required)	\$87.43
<input type="checkbox"/> PNG30306 Velcro® Positioning Belt-38"-67" (back required)	\$87.43
<input type="checkbox"/> PNG30028 Positioning Belt with Airline Style Buckle-26"-55" (back required)	\$213.86

Backs

<input type="checkbox"/> PNG50066 Flat Back-15"H	\$926.71
<input type="checkbox"/> PNG50067 Flat Back-19"H	\$977.82
<input type="checkbox"/> PNG50062 Contoured Back-15"H	\$1,019.51
<input type="checkbox"/> PNG50068 Contoured Back-19"H	\$1,070.62
<input type="checkbox"/> PNG50171 Removable Flat Back-15"H	\$1,062.55
<input type="checkbox"/> PNG50173 Removable Flat Back-19"H	\$1,113.66
<input type="checkbox"/> PNG50172 Removable Contoured Back-15"H	\$1,156.70
<input type="checkbox"/> PNG50174 Removable Contoured Back-19"H	\$1,209.16

Back Accessories

<input type="checkbox"/> P600002 Back Angle Adjustment Wheel (back required)	\$145.26
<input type="checkbox"/> P80246 Angle Locator	\$53.80

Chest Straps-back required

<input type="checkbox"/> PNG50192 Chest Strap-30"-46" (back required)	\$213.86
<input type="checkbox"/> PNG50193 Chest Strap-34"-54" (back required)	\$213.86
<input type="checkbox"/> PNG50504 Chest Strap-46"-66" (back required)	\$213.86

Upper Body Support-back required

ATTENTION - ALL OPTIONS BELOW REQUIRE THE ACCESSORIES MOUNTING BRACKET

- PNG50034 Accessories Mounting Bracket \$95.50

Chest Vests-back and accessories mounting bracket required

- PNG11212 High Mount Chest Vest Bracket \$108.95
- PNG30170 X-Style Chest Vest-11"Lx9.5"W \$399.47
- PNG30056 X-Style Chest Vest-14.5"Lx10.5"W \$399.47
- PNG30171 Y-Style Chest Vest-17.5"Lx11.5"W \$399.47
- PNG30021 Y-Style Chest Vest-20.5"Lx12.5"W \$399.47
- PNG30318 Non-Stretch X-Style Chest Vest-11"Lx9.5"W \$464.03
- PNG30322 Non-Stretch X-Style Chest Vest-14.5"Lx10.5"W \$464.03
- PNG30320 Non-Stretch Y-Style Chest Vest-17.5"Lx11.5"W \$464.03
- PNG30324 Non-Stretch Y-Style Chest Vest-20.5"Lx12.5"W \$464.03

Lateral Supports-back and accessories mounting bracket required

- PNG50535 Lateral Supports with Flat Pads-8"-16.5"W (pair) \$505.72
- PNG50536 Lateral Supports with Flat Pads-11"-19.5"W (pair) \$505.72
- PNG50224 Lateral Supports with Curved Pads-8"-16.5"W (pair) \$505.72
- PNG50176 Lateral Supports with Curved Pads-11"-19.5"W (pair) \$505.72

Head Supports-back and accessories mounting bracket required

- PNG50240 Head Support-Short \$521.86
- PNG50037 Head Support-Tall \$521.86

Other Accessories

- PNG50040 Push Handles (back required, pair) \$213.86
- P80809 Tool Pouch \$36.32

Notes:

(A) Seat Depth Range: 24"-29"
(center of knee pad to flat back)
 Measure the distance between the back of the buttocks and the back of the knee.

(B) Seat to Foot Plate Range: 15.5"-18.5"
(from seat pivot)
 Measure the distance between the foot/shoe and the knee bend.



AltimateMedical

Ph: 800.342.8968 or 507.697.6393
 Fax: 877.342.8968 or 507.697.6900
 info@easystand.com easystand.com