

Bill To <small>(Altimate Medical Supplier Only)</small>	Ship To
Purchase Order #: _____ Quote #: _____	Name: _____
Account #: _____ Ordered By: _____	Address: _____
Supplier Company: _____	Address: _____
Address: _____ City, State, Zip: _____	City, State: _____
Phone: _____ Fax: _____	Zip: _____
Email: _____ Payer Source: _____	Phone: _____
Name of ATP (Required): _____	Tag For: _____

Glider Medium User Specifications	
Approximate Height Range	4'0"-5'6"
Weight Limit	200 lbs.
Seat Depth Range (center of knee pad to flat back)	20"-25"
Seat to Foot Plate Range (from seat pivot)	13.5"-16.5"

Base

<p>EasyStand® Glider Medium \$7,246</p> <p>Price Includes:</p> <p>PNG50378 Glider Base \$3,313</p> <p>PNG50187 Glider Configuration \$3,933</p>
--

The base unit includes: Full Range Seat, Adjustable Resistance Cylinders, Positioning Bar, Two 5" Front Wheels, Two 5" Rear Locking Casters, Manual Hydraulic Actuator with Handle, Tray & Chest Pad, Glide Handles, Flip-Up Knee Pads, Adjustable Foot Plates, Transfer Handles, and Black Upholstery.

Note: The default items will be added if no selection is made.

Default items include: 1-Piece Knee Pad.

Choose actuator handle location: Right Side **default** Left Side

Seat Adjustment-CANNOT be ordered off-unit

PNG50499 Easy-Adjust Seat Depth (roller bearing allows for smoother seat depth adjustment) \$108

Lift Mechanisms

- PNG50010 Pow'r Up Lift (replaces hydraulic actuator, battery powered) \$1,222
- PNG50328 Adjustable/Removable Actuator Handle (not available with pow'r up) \$228
- PNG50044 Quad Grip Handle Extension (not available with pow'r up) \$153
- PNG50070 T-Style Handle Extension (not available with pow'r up) \$153

Foot Plate Accessories

- PNG30030 Foot Straps-15"L (length over top of foot, pair) \$46
- PNG50505 Foot Straps-21"L (length over top of foot, pair) \$46
- PNG30031 Secure Foot Straps-15"L (length over top of foot, two pair) \$117
- PNG50506 Secure Foot Straps-21"L (length over top of foot, two pair) \$117

Seats

PNG50351 Seat with Roho® Insert \$609

Seat Accessories

- PNG50345 Hip Supports (9"-15" range, pair) \$363
- PNG50374 Hip Supports (11"-17" range, pair) \$363

Handle Accessories

PNG50251 Glide Handle Extensions (pair) \$153

Tray Accessories

PNG50295 No-Tray Chest Pad \$253

Backs

- PNG50066 Flat Back-15"H \$665
- PNG50067 Flat Back-19"H \$702
- PNG50062 Contoured Back-15"H \$732
- PNG50068 Contoured Back-19"H \$769
- PNG50171 Removable Flat Back-15"H \$763
- PNG50173 Removable Flat Back-19"H \$800
- PNG50172 Removable Contoured Back-15"H \$830
- PNG50174 Removable Contoured Back-19"H \$868

Positioning and Upper Body Support

ATTENTION - All options below require a back

<input type="checkbox"/> PNG30029 Velcro® Positioning Belt-26"-55"	\$57
<input type="checkbox"/> PNG30306 Velcro® Positioning Belt-38"-67"	\$57
<input type="checkbox"/> PNG30028 Positioning Belt with Airline Style Buckle-26"-55"	\$153
<input type="checkbox"/> PNG50192 Chest Strap-30"-46"	\$142
<input type="checkbox"/> PNG50193 Chest Strap-34"-54"	\$142
<input type="checkbox"/> PNG50504 Chest Strap-46"-66"	\$142
<input type="checkbox"/> P600002 Back Angle Adjustment Wheel	\$104

ATTENTION - All options below require the accessories mounting bracket

<input type="checkbox"/> PNG50034 Accessories Mounting Bracket	\$68
<input type="checkbox"/> PNG30170 X-Style Chest Vest-11"Lx9.5"W	\$277
<input type="checkbox"/> PNG30056 X-Style Chest Vest-14.5"Lx10.5"W	\$277
<input type="checkbox"/> PNG30171 Y-Style Chest Vest-17.5"Lx11.5"W	\$277
<input type="checkbox"/> PNG30021 Y-Style Chest Vest-20.5"Lx12.5"W	\$277
<input type="checkbox"/> PNG30318 Non-Stretch X-Style Chest Vest-11"Lx9.5"W	\$323
<input type="checkbox"/> PNG30322 Non-Stretch X-Style Chest Vest-14.5"Lx10.5"W	\$323
<input type="checkbox"/> PNG30320 Non-Stretch Y-Style Chest Vest-17.5"Lx11.5"W	\$323
<input type="checkbox"/> PNG30324 Non-Stretch Y-Style Chest Vest-20.5"Lx12.5"W	\$323
<input type="checkbox"/> PNG11212 High Mount Chest Vest Bracket	\$78
<input type="checkbox"/> PNG50535 Lateral Supports with Flat Pads-8"-16.5"W (pair)	\$363
<input type="checkbox"/> PNG50536 Lateral Supports with Flat Pads-11"-19.5"W (pair)	\$363
<input type="checkbox"/> PNG50224 Lateral Supports with Curved Pads-8"-16.5"W (pair)	\$363
<input type="checkbox"/> PNG50176 Lateral Supports with Curved Pads-11"-19.5"W (pair)	\$363

Head Support

<input type="checkbox"/> PNG50240 Head Support-Short (back and accessories mounting bracket required)	\$363
<input type="checkbox"/> PNG50037 Head Support-Tall (back accessories mounting bracket required)	\$363

Other Accessories

<input type="checkbox"/> P80246 Angle Locator	\$36
<input type="checkbox"/> PNG50040 Push Handles (back and accessories mounting bracket required, pair)	\$153
<input type="checkbox"/> P80809 Tool Pouch	\$25

Notes:

AltimateMedical

Ph: 800.342.8968 or 507.697.6393

Fax: 877.342.8968 or 507.697.6900

info@easystand.com easystand.com

Email form to webquotes@easystand.com